

BOSTON OVERLAP BRACE

A Guide for Patients


 oapl.™



INTRODUCTION

Your Boston Overlap Brace "BOB" or Anterior Overlap Brace has been fitted to control the alignment and reduce the movement of your lower thoracic and lumbar spine.

Make sure you have discussed the use of your brace with both your orthotist and your doctor prior to discharge.



Consulting Orthotist

Clinic Attended

Date Attended

APPLICATION

Always wear a firm fitting shirt or singlet under the brace. Preferably one made of natural fibres such as cotton.

This brace can be applied either in sitting or lying positions. Speak to your orthotist about which is most suitable.

Lying Position

1. Log roll onto your side and have someone slide the brace under you.



2. Then roll back into the brace.

3. Locate pelvic roll and tighten all straps starting with the top one before sitting up.



Sitting Position

1. Log roll and then sit up on the edge of the bed keeping your back straight and being careful not to twist (use the technique shown to you by your physiotherapist).

2. Locate the pelvic roll in the middle of the brace on your waist (top of your pelvis) with the opening to the front.

3. Ask your orthotist to mark the top and bottom of the brace for you. Fasten the 3 straps firmly beginning with the middle strap. Breathe in and re-tension all straps.



GENERAL CARE

- The brace is to be worn at all times when you are out of bed.
- You may or may not be able to remove the brace for showering.
- Discuss this with your orthotist and doctor to clarify which you should do.
- Avoid the use of creams and powders under your brace as this may lead to skin irritations.
- Your brace may be cleaned with warm soapy water.
- Fluff should be cleaned from the Velcro straps as required to ensure that the grip remains adequate. A toothbrush can be used to do this.

Do Not

1. Heat brace
2. Grind, cut, file brace or straps. If it is uncomfortable please contact your orthotist.
3. Do not place anything under your brace as it may create pressure areas.

COMMON PROBLEMS

Heat Retention and Swelling

Unfortunately there is little you can do to combat these symptoms. Ventilation holes make little difference and often lead to pressure problems. Regular changes of the shirt/singlet worn under the brace and periods of lying down out of the brace are the best remedies.

Brace Riding Up When Sitting

Check you are placing the brace low enough when putting it on and tighten the top strap first. If it persists contact your orthotist.

Discomfort from the Pelvic Roll

It is quite common to get some discomfort in this region when first wearing your brace, especially if you have had a bone graft taken from your pelvis. Generally no modifications are required and you will become accustomed to the pressure.

Pressure Over the Top of the Thighs When Sitting

Some pressure is required in this area especially if your spinal problem is in the lower regions such as L4 or L5 levels. However, if you feel it is too much, contact your orthotist.

Pressure Under the Breasts

Ladies should always wear a bra with their brace, which should help to prevent any rubbing or pressure. If problems persist contact your orthotist so that the brace may be modified.

Straps Too Long or Short

If you gain or loose weight the straps on your brace may need some adjustment. Contact your orthotist to have this adjusted.



CONTACT DETAILS

Fitzroy Clinic

93 - 97 Webb Street
Fitzroy VIC 3065
T: 03 9419 2499
F: 03 9416 3534
mainfacility@oapl.com.au

Box Hill Clinic

Epworth Eastern
Suite 5A, Level 2, Arnold Street
Box Hill VIC 3128
T: 03 9897 1374
F: 03 9783 6944

Bundoora Clinic

Northpark Private Hospital
Greenhills Road
Bundoora VIC 3083
T: 03 9467 1222

Canberra Clinic

Mobility Matters
35 Townsville Street
Fyshwick ACT 2609
T: 02 6280 7244
F: 02 6239 1281
oapl@mobilitymatters.com.au

Clayton Clinic

281 Clayton Road
Clayton VIC 3168
T: 03 9562 9422
F: 03 9562 9605
clayton@oapl.com.au

Footscray Clinic

55 Pickett Street
Footscray VIC 3011
T: 03 9687 4589
F: 03 9687 3961
footscray@oapl.com.au

Frankston Clinic

346 Nepean Highway
Frankston VIC 3199
T: 03 9783 3866
F: 03 9783 6944
frankston@oapl.com.au

Malvern Clinic

1489 Malvern Road,
Malvern VIC 3144
T: 03 9824 5874
F: 03 9824 8205
malvernclinic@oapl.com.au

Richmond Clinic

Epworth Centre
Level 7, Suite 5, 32 Erin Street
Richmond VIC 3121
T: 03 9421 6226
F: 03 9426 4321

Ringwood Clinic

86 Mt Dandenong Road
Ringwood VIC 3135
T: 03 9879 4299
F: 03 9870 9411
ringwood@oapl.com.au

Bendigo Clinic

Anne Caudle Centre
Barnard Street
Bendigo VIC 3550
T: 03 5454 8752
F: 03 5454 8756
oapl@bendigohealth.org.au

Main Facility and Administration

93-97 Webb Street

Fitzroy Victoria 3065

T: 03 9419 2499

F: 03 9416 3543

mainfacility@oapl.com.au

ACN 006 520 314

ABN 30 006 520 314

www.oapl.com.au

Version 1.09

OAPT0017

