

Juzo® Arion Slip-on and Taking-off Aids

For Juzo® Compression Garments

The expanded range of Juzo Arion Products make it simple and easy to put on and take off medical compression garments such as stockings, pantyhose, or sleeves. Thanks to the very smooth material, putting on and taking off a compression garment requires very little physical effort on the part of the patient or nursing staff. This increases patient compliance when a compression garment needs to be worn on a regular basis.

Characteristics

- very smooth, durable and tear-resistant material
- foldable and therefore easy to carry
- does not damage fragile skin
- can be used for scar therapy and lymphology
- washable (hand wash)

	Compression sleeves		Compression stockings and pantyhose		Above-knee prosthesis
		suitable for children	open toe	closed toe	with suction socket
Slip-on and taking-off aids					
Juzo Arion Sim-Slide			X		
Juzo Arion Easy-Slide Leg + Easy-Off			X	X**	
Slip-on aids					
Juzo Arion Easy-Slide Arm	X	X*			
Juzo Arion Easy-Slide Leg			X		
Juzo Arion Magnide				X	
Juzo Arion SlideX			X	X	
Juzo Arion Proth-Pro					X

* Size XS also suitable as a slip-on aid leg for children.

** Only possible with the Juzo Arion Easy-Off.

**International
Julius Zorn GmbH**
Juliusplatz 1
86551 Aichach
Germany
info@juzo.com

www.juzo.com

Juzo UK Ltd
South Court
Sharston Rd
Sharston
Manchester M22 4 SN
Great Britain
info@juzo.co.uk

www.juzo.co.uk



ARION SLIP-ON AND TAKING-OFF AIDS

For Juzo® Compression Garments

Juzo® Arion Easy-Slide Arm

Very easy slipping-on of compression sleeves in only 5 steps



Slip your arm into the Juzo Arion Easy-Slide Arm.

Now pull the compression sleeve over the slip-on aid on your arm and put the silicone border at the right place.

Finally hang the loop at the end of the Juzo Arion Easy-Slide Arm onto e.g. a door handle or hold it with the other hand. Pull the slip-on aid slowly out from under the compression sleeve. Finished!



EN • 14004803/3/3 • W362 • 05/2016

Juzo® Arion Slip-on and Taking-off Aids

NEW

Juzo® Arion Sim-Slide

Leg slip-on and taking-off aid for compression stockings/pantyhose with **open** toe

- smooth, two-layer material
- increases the independence and therapy compliance of patients
- reduces the physical exertion required by the patient or nursing staff
- simply fasten to the ankle using the elastic strap



NEW

Juzo® Arion Easy-Slide Leg + Easy-Off

Leg slip-on and taking-off aid for compression stockings/pantyhose with **open** toe

- smooth, two-layer material
- increases the independence and therapy compliance of patients
- reduces the physical exertion required by the patient or nursing staff
- taking-off aid is also suitable for closed toe
- with velcro fastener for easy fixation above the ankle



Juzo® Arion Slip-on Aids

Juzo® Arion Easy-Slide Arm

Arm slip-on aid for compression sleeves

- smooth, two-layer material
- increases the independence and therapy compliance of patients
- reduces the physical exertion required by the patient or nursing staff
- for children size XS also suitable as a slip-on aid for legs



Juzo® Arion Easy-Slide Leg

Leg slip-on aid for compression stockings/pantyhose with **open** toe

- smooth, two-layer material
- increases the independence and therapy compliance of patients
- reduces the physical exertion required by the patient or nursing staff



Juzo® Arion Magnide

Leg slip-on aid for compression stockings/pantyhose with **closed** toe

- smooth material
- increases the independence and therapy compliance of patients
- reduces the physical exertion required by the patient or nursing staff
- with innovative magnetic closure



Juzo® Arion SlideX

Extension for Juzo Arion Easy-Slide Leg and Juzo Arion Magnide

- smooth material
- reduces the physical exertion required by the patient or nursing staff
- especially for large leg circumferences
- can only be used in combination with Juzo Arion Easy-Slide Leg and Juzo Arion Magnide



NEW

Juzo® Arion Proth-Pro

Slip-on aid for thigh prosthesis with suction socket

- smooth material
- reduces the physical exertion required by the patient
- with extender for even easier use



Juzo® Arion Sim-Slide

Art. no.: 1705

Size	Shoe size
M	38 - 41
L	42 - 45
XL	≥ 46

Juzo® Arion Easy-Slide Arm

Art. no.: 1701

Size	Arm length
S	39.5 cm
M	55.5 cm
L	62.5 cm

Juzo® Arion Easy-Slide Leg

Art. no.: 1702

Size	Shoe size
S	34 - 37
M	38 - 41
L	42 - 45
XL	≥ 46

Juzo® Arion SlideX

Art. no.: 1704

Size	Thigh circumference
M	≤ 40 cm
L	41 - 60 cm
XL	61 - 100 cm

Juzo® Arion Easy-Slide Leg + Easy-Off

Art. no.: 1707

Size	Shoe size
S	34 - 37
M	38 - 41
L	42 - 45
XL	≥ 46

Size	Arm / leg length
XS	25.5 cm

FOR CHILDREN

Juzo® Arion Magnide

Art. no.: 1703

Size	Shoe size
M	≤ 40
L	41 - 45
XL	≥ 46

The size XL is also suitable for wide feet with a circumference of up to 28 cm.

Juzo® Arion Proth-Pro

Art. no.: 1706

Size	Circumference	Length
XS	≤ 40 cm	12 cm
S	41 - 45 cm	13 cm
M	46 - 58 cm	13 cm
L	59 - 70 cm	15.5 cm
XL	71 - 80 cm	16 cm
XXL	≥ 81 cm	17 cm